



PRODUCT CODE: 281380

Acidimetric Liquor titrated for volumetric analysis

for determination of oils and fats acidity. 1 ml corresponds to 0,028245 g of oleic acid

TARIC 3822 90 00 00

PHYSICAL DATA: liquid, Clear, Colourless, Miscible with water • D 20/4 1,002 •

BIBLIOGRAPHY:

HAZARDOUS: LD50 oral rat 273 mg/kg •



H: H319 • H315 •

P: P264 • P280 • P302+P352 • P305+P351+P338 • P321 • P501 • P332+P313 • P337+P313 • P362 •

TRANSPORT REGULATIONS: UN: 1719 • ADR: 8/III • IMDG: 8/III • IATA: 8/III • PAX: 852 • CAO: 856 • (E) •

WEIGHT/VOLUME INFORMATION: 1l-1,002 kg 1kg-0,998 l

COMPOSITION:

Potassium Hydroxide 50% w/w
Water (s.q.m)

1,4 g
100 ml