



PRODUCT CODE: 281381

Acidimetric Liquor titrated for volumetric analysis

for determination of oils and fats acidity in commercial grades. For 10 g of oil, 1 ml corresponds to 1° of acidity (=0,1 g of oleic acid)

TARIC 3822 90 00 00

PHYSICAL DATA: liquid, Clear, Colourless, Miscible with water • D 20/4 1,01 •

BIBLIOGRAPHY:

HAZARDOUS: LD L0 oral rbt 500 mg/kg • LD50 ipr mus 40 mg/kg TLV-TWA (NaOH) 2 mg/m³ VLA-EC (NaOH) 2 mg/m³



H: H319 • H315 • H290 •

P: P264 • P280 • P302+P352 • P305+P351+P338 • P321 • P501 • P332+P313 • P337+P313 • P362 •

TRANSPORT REGULATIONS: UN: 1719 • ADR: 8/III • IMDG: 8/III • IATA: 8/III • PAX: 852 • CAO: 856 • (E) •

WEIGHT/VOLUME INFORMATION: 1l~1,010 kg 1kg~0,99 l

COMPOSITION:

Sodium Hydroxide 50% sol. w/w	2,06 ml
Water (s.q.m)	100 ml