



PRODUCT CODE: 281572

Hanus Reagent 0.1 mol/l (0.2N) for volumetric analysis
for determination of iodine index Keep over 18°C.

TARIC 3822 90 00 00

PHYSICAL DATA: liquid, reddish yellow, Miscible with water • D 20/4 1,069 • pH1 • Flash P.: ~40 °C •

BIBLIOGRAPHY: AOAC , 41 , 6 (1995), Métodos Oficiales de Análisis I , 72 (1993) •

HAZARDOUS: CEE (referring to acetic acid). 10 ml/m3 25 mg/m3 VLA-EC (referring to acetic acid). 15 ppm 37 mg/m3 VLA-ED (referring to acetic acid). 10 ppm 25 mg/m3



H: H226 • H314 • H290 • H312 •

P: P280 • P303+P361+P353 • P305+P351+P338 • P310 • P321 • P501 • P362+P364 • P405 •

TRANSPORT REGULATIONS: UN: 2789 • ADR: 8(3)/II • IMDG: 8(3)/II • IATA: 8(3)/II • PAX: 851 • CAO: 855 • (E) •

WEIGHT/VOLUME INFORMATION: 1l-1,069 kg 1kg-0,935 l

COMPOSITION:

Iodine	1,32 g
Bromine	0,812 g
Acetic Acid glacial s.q.m	100 ml