



PRODUCT CODE: 281730

Buffer Solution pH 10 for volumetric analysis
for complexometry

TARIC 3822 90 00 00

PHYSICAL DATA: liquid, Clear, Colourless, Miscible with water • D 20/4 0,977 • pH10 •

BIBLIOGRAPHY:

HAZARDOUS: / Data not available. • VLA-EC (NH3) 25 mg/m3 VLA-EC (NH3) 35 ppm VLA-ED (NH3) 25 ppm18 mg/m3



H: H314 • H411 •

P: P260 • P264 • P273 • P280 • P301+P330+P331 • P501 • P303+P361+P353 • P304+P340 • P305+P351+P338 • P310 • P321 • P338 • P363 • P391 • P405 •

TRANSPORT REGULATIONS: UN: 1719 • ADR: 8/III • IMDG: 8/III • IATA: 8/III • PAX: 852 • CAO: 856 • (E) •

WEIGHT/VOLUME INFORMATION: 1l-0,977 kg 1kg-1,024 l

COMPOSITION:

Ammonium Chloride	6,75 g
Ammonia 30%	35 ml
Water (s.q.m)	100 ml

SPECIFICATIONS:

Identity :	
Identity	passes test
pH	10,0 - 10,5

