



PRODUCT CODE: 282430

Indicator Tashiro 4.4. Mixed (Methyl Red-Methylene Blue) for volumetric analysis
for ammonia titrations. pH 4,4 red violet; 5,8 green

TARIC 3822 90 00 00

PHYSICAL DATA: liquid, Clear, grey to brownish, Miscible with alcohol and water • D 20/4 0,828 • Flash P.: ~14 °C
• Ign. T.: ~425 °C •

BIBLIOGRAPHY:

HAZARDOUS: LD50 oral rat 7.060 mg/kg • VLA-ED (ethanol) 1.000 ppm 1.910 mg/m³



H: H225 • H319 •

P: P210 • P233 • P240 • P241 • P242 • P501 • P243 • P280 • P303+P361+P353 • P370+P378 • P403+P235 •

TRANSPORT REGULATIONS: UN: 1993 • ADR: 3/II • IMDG: 3/II • IATA: 3/II • PAX: 353 • CAO: 364 • (D/E) •

WEIGHT/VOLUME INFORMATION: 1l~0,828 kg 1kg~1,208 l

COMPOSITION:

Methyl Red	100 mg
Methylene Blue	50 mg
Water	10,4 ml
Ethanol absolute s.q.m	100 ml

SPECIFICATIONS:

Identity :	
Identity	passes test
pH transition interval:	
violet red	4,4

green

5,8