



PRODUCT CODE: 285483

Ethanol / Diethyl Ether 1:1 v/v + Phenolphthalein 0,1 % for acidity determination in olive oil

TARIC 3822 90 00 00

PHYSICAL DATA: liquid, Colourless, Miscible with alcohol • D 0,765 • Flash P.: <10,0 °C(closed cup) •

BIBLIOGRAPHY: Métodos Analíticos en Alimentaria (Métodos Oficiales de Análisis) Aceites y grasas, Panreac, Métodos Oficiales de Análisis de Aceites y Grasas I , 51 (1993), Ministerio de Agricultura, Pesca y Alimentación, Standard Methods for the Analysis of Oils, Fats and Soaps, II.D.I (1964) , IUPAC •

HAZARDOUS: LD L0 oral man 260 mg/kg • LD50 oral rat 1.215 mg/kg • LC50 rat 73000ppm (V) / 2h • VLA-EC (diethyl ether) 550 ppm1.540 mg/m3 VLA-ED (ethanol) 1.000 ppm1.910 mg/m3 VLA-ED (diethyl ether) 400 ppm1.230 mg/m3



H: H224 • EUH019 • H319 • H336 •

P: P210 • P280 • P303+P361+P353 • P305+P351+P338 • P405 • P501 •

TRANSPORT REGULATIONS: UN: 1993 • ADR: 3/I • IMDG: 3/I • IATA: 3/I • PAX: 351 • CAO: 361 • (D/E) •

WEIGHT/VOLUME INFORMATION: 1l-0,765 kg 1kg-1,307 l

COMPOSITION:

Phenolphthalein	100 mg
Diethyl Ether	50 ml
Ethanol	50 ml

SPECIFICATIONS:

Identity :	
Identity	passes test
Density at 20/4	0,755 - 0,762

