



PRODUCT CODE: 285813

AQUAMETRIC Composite 2 for Karl Fischer volumetric analysis

for volumetric Karl Fischer titrations. One-component reagent. 1 ml corresponds to 2-2,4 mg of H₂O

TARIC 3822 90 00 00

PHYSICAL DATA: liquid, D 20/4 1,11 • B.P.: 194 °C • Flash P.:90 °C •

BIBLIOGRAPHY:

HAZARDOUS:



H: H314 • H372 • H360D • H302 • H332 •

P: P260 • P261 • P264 • P270 • P271 • P273 • P280 • P301+P312 • P301+P330+P331 • P302+P352 • P303+P361+P353 • P304+P340 • P305+P351+P338 • P391 • P501 •

TRANSPORT REGULATIONS: UN: 1760 • ADR: 8/III • IMDG: 8/III • IATA: 8/III • PAX: 852 • CAO: 856 • (E) •

WEIGHT/VOLUME INFORMATION: 1l-1,110 kg 1kg-0,901 l

SPECIFICATIONS:

1 ml corresponds to:

2 - 2,4 mg of H₂O*

*** At the moment of the batch analysis.**

Protect from moisture.

Ed.: 5 . Vig.: 19.02.2019 .

Prod.: 285813