



PRODUCT CODE: 285814

AQUAMETRIC Composite 5K for Karl Fischer volumetric analysis

for volumetric Karl Fischer titrations in ketones and aldehydes. One-component reagent. 1 ml corresponds to 5 - 5,5 mg of H₂O

TARIC 3822 90 00 00

PHYSICAL DATA: liquid, Clear, Dark brown, D 20/4 1,17 • B.P.: 194 °C • Flash P.:90 °C •

BIBLIOGRAPHY:

HAZARDOUS:



H: H314 • H302+H332 • H360D • H372 •

P: P280 • P303+P361+P353 • P305+P351+P338 • P310 • P501 • P405 • P321 •

TRANSPORT REGULATIONS: UN: 1760 • ADR: 8/III • IMDG: 8/III • IATA: 8/III • PAX: 852 • CAO: 856 • (E) •

WEIGHT/VOLUME INFORMATION: 1l~1,170 kg 1kg~0,855 l

SPECIFICATIONS:

1 ml corresponds to:

5 - 5,5 mg of H₂O (20°C)*

*** At the moment of the batch analysis.**

Protect from moisture.

