



PRODUCT CODE: 285817

AQUAMETRIC Solvent for Karl Fischer volumetric analysis
for volumetric Karl Fischer titrations. Use with AQUAMETRIC Titrant

TARIC 3822 90 00 00

PHYSICAL DATA: liquid, Colourless, D 20/4 0,891 • B.P.: 65 °C • Flash P.:11 °C •

BIBLIOGRAPHY:

HAZARDOUS: VLA-EC (methanol) 250 ppm333 mg/m3 VLA-ED (methanol) 200 ppm266 mg/m3



H: H225 • H331 • H311 • H301 • H314 • H370 •

P: P210 • P233 • P240 • P241 • P242 • P501 • P243 • P260 • P261 • P264 • P270 • P271 • P280 • P301+P310 • P301+P330+P331 • P302+P352 • P303+P361+P353 • P304+P340 • P305+P351+P338 • P307+P311 • P310 • P311 • P312 • P321 • P322 • P330 • P338 • P361 • P363 • P370+P378 • P403+P233 • P403+P235 • P405 •

TRANSPORT REGULATIONS: UN: 3286 • ADR: 3(6.1)(8)/II • IMDG: 3(6.1)(8)/II • IATA: 3(6.1)(8)/II • PAX: 352 • CAO: 363 • (D/E) •

WEIGHT/VOLUME INFORMATION: 1l~0,877 kg 1kg~1,14 l

SPECIFICATIONS:

Suitability: for H₂O determination

passes test

Close the bottle immediately after use.

Protect from moisture.

Ed.: 4 . Vig.: 21.11.2018 .

Prod.: 285817