



PRODUCT CODE: 285819

AQUAMETRIC Solvent CM for Karl Fischer volumetric analysis
for volumetric Karl Fischer titrations in oils and fats. Use with AQUAMETRIC Titrant

TARIC 3822 90 00 00

PHYSICAL DATA: liquid, Colourless, D 20/4 1,277 • B.P.: 60 - 65 °C • Flash P.: >23 °C •

BIBLIOGRAPHY:

HAZARDOUS: VLA-ED (CHCl₃) 10 ppm/50 mg/m³



H: H226 • H302+H312 • H331 • H315 • H318 • H372 • H370 • H360D • H351 •

P: P280 • P303+P361+P353 • P305+P351+P338 • P310 • P321 • P330 • P501 • P405 • P362+P364 •

TRANSPORT REGULATIONS: UN: 1992 • ADR: 3(6.1)/III • IMDG: 3(6.1)/III • IATA: 3(6.1)/III • PAX: 352 • CAO: 364 • (D/E) •

WEIGHT/VOLUME INFORMATION: 1l-1,277 kg 1kg-0,783 l

SPECIFICATIONS:

Suitability: for H₂O determination
Close the bottle immediately after use.

passes test

Protect from moisture.

