



PRODUCT CODE: 285821

AQUAMETRIC Working Medium for Karl Fischer volumetric analysis

for volumetric Karl Fischer titrations in ketones and aldehydes. Use with AQUAMETRIC Composite 5K

TARIC 3822 90 00 00

PHYSICAL DATA: liquid, Clear, Colourless, D 20/4 1,356 • B.P.: 60 °C •

BIBLIOGRAPHY:

HAZARDOUS: LD50 oral rat 76 mg/kg • LC50 rat 0.11mg/l / 4h • LD50 skn rbt 68 mg/kg •



H: H330 • H311 • H301 • H315 • H351 • H372 • H319 • H361d •

P: P201 • P202 • P260 • P262 • P264 • P501 • P270 • P271 • P280 • P281 • P284 • P301+P310 • P302+P350 • P302+P352 • P304+P340 • P308+P313 • P310 • P314 • P320 • P321 • P322 • P330 • P332+P313 • P361 • P362 • P363 • P403+P233 • P405 •

TRANSPORT REGULATIONS: UN: 2810 • ADR: 6.1/II • IMDG: 6.1/II • IATA: 6.1/II • PAX: 654 • CAO: 662 • (D/E) •

WEIGHT/VOLUME INFORMATION: 1l-1,356 kg 1kg-0,737 l

SPECIFICATIONS:

Suitability: for H₂O determination

passes test

Protect from moisture.

