



PRODUCT CODE: 286180

AQUAMETRIC Coulomat AG for Karl Fischer coulometric analysis

anolyte for the coulometric determination of water. For cells with or without diaphragm.

TARIC 3822 90 00 00

PHYSICAL DATA: liquid, yellow to brown, Very slightly soluble in water. D 20/4 0,915 • B.P.: 64 - 65 °C • Flash P.:12 °C •

BIBLIOGRAPHY:

HAZARDOUS:



H: H225 • H331 • H311 • H301 • H314 • H370 • H372 • H360D •

P: P201 • P202 • P210 • P233 • P240 • P501 • P241 • P242 • P243 • P260 • P261 • P264 • P270 • P271 • P280 • P281 • P301+P310 • P301+P330+P331 • P302+P352 • P303+P361+P353 • P304+P340 • P305+P351+P338 • P307+P311 • P308+P313 • P310 • P311 • P312 • P321 • P322 • P330 • P338 • P361 • P363 • P370+P378 • P403+P233 • P403+P235 • P405 •

TRANSPORT REGULATIONS: UN: 3286 • ADR: 3(6.1)(8)/II • IMDG: 3(6.1)(8)/II • IATA: 3(6.1)(8)/II • PAX: 352 • CAO: 363 • (D/E) •

WEIGHT/VOLUME INFORMATION: 1l~0,915 kg 1kg~1,093 l

SPECIFICATIONS:

Suitability: for H₂O determination

passes test

Contains: Methanol, Imidazole, Sulfur Dioxide, Iodine, Diethanolamine.

