



PRODUCT CODE: 287192

AQUAMETRIC Coulomat CG for Karl Fischer coulometric analysis

Catholyte for the coulometric determination of water. For cells with diaphragm. Use with AQUAMETRIC Coulomat A or Coulomat AG

TARIC 3822 90 00 00

PHYSICAL DATA: Very slightly soluble in water. D 20/4 0,915 •

BIBLIOGRAPHY:

HAZARDOUS:



H: H225 • H331 • H311 • H301 • H314 • H370 • H360D • H373 •
P: P201 • P202 • P210 • P233 • P240 • P501 • P241 • P242 • P243 • P260 • P261 • P264 • P270 • P271 • P280 • P281 •
P301+P310 • P301+P330+P331 • P302+P352 • P303+P361+P353 • P304+P340 • P305+P351+P338 • P307+P311 •
P308+P313 • P310 • P311 • P312 • P321 • P322 • P330 • P338 • P361 • P363 • P370+P378 • P403+P233 • P403+P235 • P405 •

TRANSPORT REGULATIONS: UN: 3286 • ADR: 3(6.1)(8)/II • IMDG: 3(6.1)(8)/II • IATA: 3(6.1)(8)/II • PAX: 352 •
CAO: 363 • (D/E) •

WEIGHT/VOLUME INFORMATION: 1l~0,915 kg 1kg~1,092 l

SPECIFICATIONS:

Suitability: for H₂O determination passes test
Contains: Methanol, Imidazole, Sulfur Dioxide, Iodine,
Diethanolamine.

