



PRODUCT CODE: 303114

Sulfuric Acid 0.05 mol (4.904g H₂SO₄) to prepare 1l of 0.1N solution

H₂O₄S

H₂SO₄

M.= 98,08

CAS [7664-93-9]

EINECS 231-639-5

TARIC 2807 00 00 00

PHYSICAL DATA: Plastic or glass ampoules containing clear colourless liquid, Miscible with water • D 20/4 1,07 • pH1 •

BIBLIOGRAPHY:

HAZARDOUS: C.E: 016-020-00-8 • VLA-EC (H₂SO₄) 3 mg/m³ VLA-ED (H₂SO₄) 1 mg/m³



H: H290 •

P: P234 • P390 • P406 •

TRANSPORT REGULATIONS: UN: 2796 • ADR: 8/II • IMDG: 8/II • IATA: 8/II • PAX: 851 • CAO: 855 • (E) •

OBSERVATIONS: Product controlled as a drug precursor. •

SPECIFICATIONS:

Titer at 20°C

0,997 - 1,003

Traceability

See certificate