



**PRODUCT CODE: 303119**

**Iodine 0.05 mol (12.690g I<sub>2</sub>) to prepare 1l of 0.1N solution**

---

CAS [7553-56-2]

EINECS 231-442-4

TARIC 2801 20 00 00

**PHYSICAL DATA:** Plastic or glass ampoules containing clear brown liquid, Miscible with water • D 1,38 • pH6,7 •

**BIBLIOGRAPHY:**

---

**HAZARDOUS:** LD L0 oral man 28 mg/kg • LD50 oral rat 14.000 mg/kg • VLA-EC (I<sub>2</sub>) 0,1 ppm1 mg/m<sup>3</sup>



H: H373 •

P: P260 • P314 • P501 •

---

**SPECIFICATIONS:**

Titer at 20°C

0,997 - 1,003

Traceability

See certificate