



**PRODUCT CODE: 303124**

**Potassium Permanganate 0.02 mol (3.161g KMnO<sub>4</sub>) to prepare 1l of 0.1N solution**

---

KMnO<sub>4</sub>

KMnO<sub>4</sub>

M.= 158,04 CAS [7722-64-7] EINECS 231-760-3

TARIC 2841 61 00 00

**PHYSICAL DATA:** Plastic or glass ampoules containing clear dark violet liquid, Miscible with water · D 1,01 · pH8 ·

**BIBLIOGRAPHY:**

**HAZARDOUS:**



H: H411 · H361d ·

P: P273 · P391 · P501 ·

**TRANSPORT REGULATIONS:** UN: 3082 · ADR: 9/III · IMDG: 9/III · IATA: 9/III · PAX: 964 · CAO: 964 · (-) ·

**OBSERVATIONS:** Product controlled as a drug precursor. ·

**SPECIFICATIONS:**

Titer at 20C  
Traceability

0,997 - 1,003  
NIST

Ed.: 5 . Vig.: 03.02.2015 .

Prod.: 303124