



PRODUCT CODE: 367170

Acetonitrile with 0.1% (v/v) of Trifluoroacetic Acid for HPLC

TARIC 3822 90 00 00

PHYSICAL DATA: liquid, Flash P.:2 °C •

BIBLIOGRAPHY:

HAZARDOUS:



H: H225 • H302 • H332 • H312 • H319 •

P: P210 • P233 • P240 • P241 • P242 • P501 • P280 • P303+P361+P353 • P370+P380 • P403+P235 • P243 • P261 • P264 • P270 • P271 • P301+P312 • P302+P352 • P304+P340 • P305+P351+P338 • P312 • P330 • P322 • P337+P313 • P363 •

TRANSPORT REGULATIONS: UN: 1648 • ADR: 3/II • IMDG: 3/II • IATA: 3/II • PAX: 353 • CAO: 364 • (D/E) •

WEIGHT/VOLUME INFORMATION: 1l-0,782 kg 1kg-1,278 l

SPECIFICATIONS:

Assay (as Trifluoroacetic Acid)

0,095 - 0,105 % (v/v)

This product has been manufactured with Acetonitrile for UHPLC (code 221881) and Trifluoroacetic Acid for UV (code 363317)