



PRODUCT CODE: 413794

Peptone Water (Dehydrated Culture Media) for microbiology

Physical Data: Dehydrated culture media of white-cream colour · Preparation: Make a suspension of 15 g in 1 l of distilledwater. Stir until completely dissolved. Distribute into tubes or flasks and sterilize at 121°C for 15 minutes ·

TARIC 3821 00 00 00

SYNONYMS: Tryptone Broth, Tryptone Water Broth•

PHYSICAL DATA:

BIBLIOGRAPHY: M. R. Pascual Anderson (2000) Microbiología alimentaria. Metodología analítica para alimentos y bebidas, Atlas R.M. & Parks L.C., Handbook of Microbiological Media , 1843 (2004) •

HAZARDOUS:

OBSERVATIONS: Storage in a cool and dry place.

Composition (g/l):

Tryptone	10,0
Sodium Chloride	5,0
pH: 7.2±0.2	