



**PRODUCT CODE: 493795**

**Buffered Peptone Water (ISO 6579, ISO 22964, ISO 6887, DIN 10181, 10160) (Prepared Bottles) for microbiology**

Physical Data: Prepared culture media of light yellow colour · Procedure: Dilute or inoculate the sample in the broth ·

---

TARIC 3821 00 00 00

**PHYSICAL DATA:**

**BIBLIOGRAPHY:** ISO 6579, ISO 22964, ISO 6887, DIN 10181, 10160 ·

**HAZARDOUS:**

**OBSERVATIONS:** Storage in a cool and dry place.

**Composition (g/l):**

Casein Peptone	10,0
Sodium Chloride	5,0
di-Sodium Phosphate 12-hydrate	9,0
mono-Potassium Phosphate	1,5
pH: 7.0±0.2	