



PRODUCT CODE: 786329

Sulfate standard solution SO₄=1.000 g/l for IC
(Na₂SO₄ in H₂O) for ionic chromatography

TARIC 3822 90 00 00

The solution of the element at the concentration given about, is NIST standard traceable .

PHYSICAL DATA: liquid, Clear, Miscible with water • D 20/4 1 •

BIBLIOGRAPHY:

HAZARDOUS:

WEIGHT/VOLUME INFORMATION: 1l~1,00 kg 1kg~1,00 l

SPECIFICATIONS:

Composition 1,479 g Na₂SO₄ / l H₂O

Concentration (as SO₄)

Uncertainty

Traceability

Method of analysis

0,990-1,010 g/l

See certificate

NIST

IC