

Specification

Nonfat dried milk powder

A0830

Physical Description:	Solid
Product Code:	A0830
Product Name:	Nonfat dried milk powder
Specifications:	Ash: max. 8.5 % Total viable counts: max. 30000 CFU/g pH (20 %; H₂O; 20°C): 6.5 ± 0.2 Protein: 33.0 - 40.0 % Fat: max. 1.0 % Lactose: 45.0 - 56.5 %
WGK:	1
Storage:	RT
EINECS:	271-045-3
CS:	04021011 (for package sizes ≤ 2.5 kg) 04021019 (for package sizes > 2.5 kg)

AppliChem GmbH

Ottoweg 4 • D-64291 Darmstadt • Phone +49 6151 9357 0 • Fax +49 6151 9357 11 • info.de@itwreagents.com • www.itwreagents.com
 CEO Joan Roget • Commerzbank Darmstadt • Bank 508 800 50 • Account 0186989900 IBAN DE24 5088 0050 0186 9899 00 • Swiftcode DRESDEFF508 • Finanzamt Darmstadt 007 225 00192 • Register court Darmstadt HRB Nr. 7340