

Specification

PBS buffer (10X Dulbecco's) - Powder

A0965

Product Code:	A0965
Product Name:	PBS buffer (10X Dulbecco's) - Powder
Specifications:	<p>Osmolality (1 X PBS buffer): 262 - 290 mOsm/kg H₂O</p> <p>pH (1 X PBS buffer; 25°C): 7.2 - 7.8</p> <p>Composition:</p> <p>Potassium chloride: 2000.00 mg/L</p> <p>Potassium dihydrogen phosphate: 2000.00 mg/L</p> <p>Sodium chloride: 80000.00 mg/L</p> <p>di-Sodium hydrogen phosphate anhydrous: 11500.00 mg/L</p>
WGK:	1
Storage:	RT
CS:	38210000

AppliChem GmbH