

Specification

YT-Medium - Powder (2X)

A0981

Product Code:	A0981
Product Name:	YT-Medium - Powder (2X)
Short Description:	usage: to set up a 2X concentrated liquid medium
Specifications:	<p>Osmolality: 286 - 316 mOsm/kg H₂O</p> <p>pH (25°C): 6.7 - 7.3</p> <p>Composition:</p> <p>Yeast extract: 10 g/L</p> <p>Sodium chloride: 5 g/L</p> <p>Tryptone: 16 g/L</p>
WGK:	1
Storage:	RT
CS:	38210000

AppliChem GmbH