

Specification

YT-Medium - Powder (2X)

A0981

| | |
|---------------------------|---|
| Product Code: | A0981 |
| Product Name: | YT-Medium - Powder (2X) |
| Short Description: | usage: to set up a 2X concentrated liquid medium |
| Specifications: | <p>pH (25°C): 6.7 - 7.3</p> <p>Osmolality: 286 - 316 mOsm/kg H₂O</p> <p>Composition:</p> <p>Yeast extract (A1552): 10 g/L</p> <p>Sodium chloride (A1149): 5 g/L</p> <p>Tryptone (A1553): 16 g/L</p> |
| WGK: | 1 |
| Storage: | RT |
| CS: | 38210000 |
| | |

AppliChem GmbH