

Specification

Bis-Tris for buffer solutions

A1025

Physical Description:	Solid
Product Code:	A1025
Product Name:	Bis-Tris for buffer solutions
Specifications:	Assay (titr.): min. 99 % pH (1 %; H₂O): 8.8 - 9.6 Heavy metals (as Pb): max. 0.0003 % Loss on drying: max. 1 % Fe: max. 0.0005 % A (1 cm/0.1 M in H₂O) 340 nm: max. 0.02 280 nm: max. 0.04
WGK:	1
Storage:	RT
Molecular Formula:	C ₈ H ₁₉ NO ₅
M:	209.24 g/mol
CAS:	6976-37-0
EINECS:	230-237-7
CS:	29221900

AppliChem GmbH

Ottoweg 4 • D-64291 Darmstadt • Phone +49 6151 9357 0 • Fax +49 6151 9357 11 • info.de@itwreagents.com • www.itwreagents.com
 CEO Joan Roget • Commerzbank Darmstadt • Bank 508 800 50 • Account 0186989900 IBAN DE24 5088 0050 0186 9899 00 • Swiftcode
 DRESDEFF508 • Finanzamt Darmstadt 007 225 00192 • Register court Darmstadt HRB Nr. 7340