

Specification

HEPES for buffer solutions

A1069

Physical Description:	Solid
Product Code:	A1069
Product Name:	HEPES for buffer solutions
Specifications:	Assay (titr.): min. 99.5 % pH (1 %; H₂O): 4.7 - 6.0 Heavy metals (as Pb): max. 0.001 % Water: max. 0.5 % Chloride: max. 0.05 % Sulfate: max. 0.05 %
WGK:	1
Storage:	RT
Molecular Formula:	C ₈ H ₁₈ N ₂ O ₄ S
M:	238.31 g/mol
CAS:	7365-45-9
EINECS:	230-907-9
CS:	29335995

AppliChem GmbH

Ottoweg 4 • D-64291 Darmstadt • Phone +49 6151 9357 0 • Fax +49 6151 9357 11 • info.de@itwreagents.com • www.itwreagents.com
 CEO Joan Roget • Commerzbank Darmstadt • Bank 508 800 50 • Account 0186989900 IBAN DE24 5088 0050 0186 9899 00 • Swiftcode
 DRESDEFF508 • Finanzamt Darmstadt 007 225 00192 • Register court Darmstadt HRB Nr. 7340