

Specification

PIPES for buffer solutions

A1079

Physical Description:	Solid
Product Code:	A1079
Product Name:	PIPES for buffer solutions
Specifications:	Assay (titr.): min. 99 % Heavy metals (as Pb): max. 0.0005 % Water: max. 1 % A (1 cm/0.1 M in 0.1 M NaOH) 260 nm: max. 0.05 280 nm: max. 0.04
WGK:	1
Storage:	RT
Molecular Formula:	$C_8H_{18}N_2O_6S_2$
M:	302.37 g/mol
CAS:	5625-37-6
EINECS:	227-057-6
CS:	29335995

AppliChem GmbH

Ottoweg 4 • D-64291 Darmstadt • Phone +49 6151 9357 0 • Fax +49 6151 9357 11 • info.de@itwreagents.com • www.itwreagents.com
 CEO Joan Roget • Commerzbank Darmstadt • Bank 508 800 50 • Account 0186989900 IBAN DE24 5088 0050 0186 9899 00 • Swiftcode DRESDEFF508 • Finanzamt Darmstadt 007 225 00192 • Register court Darmstadt HRB Nr. 7340