

Specification

Agarose MP

A1091

Physical Description:	Solid
Product Code:	A1091
Product Name:	Agarose MP
Specifications:	DNases/RNases/Proteases: not detectable Ash: $\leq 0.25 \%$ Electroendoosmosis (EEO): ≤ 0.12 Moisture: $\leq 7 \%$ Gel point: $36 \pm 1.5^{\circ}\text{C}$ Gel strength (1 %): $\geq 1800 \text{ g/cm}^2$ Gel strength (1.5 %): $\geq 3200 \text{ g/cm}^2$ Melting point: $88 \pm 1.5^{\circ}\text{C}$ Sulfate: $\leq 0.12 \%$
WGK:	1
Storage:	RT
CAS:	9012-36-6
EINECS:	232-731-8
CS:	39139000

AppliChem GmbH

Ottoweg 4 • D-64291 Darmstadt • Phone +49 6151 9357 0 • Fax +49 6151 9357 11 • info.de@itwreagents.com • www.itwreagents.com
 CEO Joan Roget • Commerzbank Darmstadt • Bank 508 800 50 • Account 0186989900 IBAN DE24 5088 0050 0186 9899 00 • Swiftcode DRESDEFF508 • Finanzamt Darmstadt 007 225 00192 • Register court Darmstadt HRB Nr. 7340