

Specification

L-Hydroxyproline pure

A1705

Solubility:	361 g/L (H ₂ O)
Physical Description:	Solid
Product Code:	A1705
Product Name:	L-Hydroxyproline pure
Headline Comment:	• All amino acids from AppliChem are of non-animal origin!
Specifications:	Assay (titr.): 99.0 - 101.0 % α 20°C/D; 5 %, H ₂ O: -74° - -77° Heavy metals (as Pb): max. 0.001 % Sulfated ash: max. 0.1 % Loss on drying: max. 0.2 % Other amino acids: max. 0.5 % Ammonium: max. 0.02 % Chloride: max. 0.02 % Sulfate: max. 0.02 %
WGK:	1
Storage:	RT
Molecular Formula:	C ₅ H ₉ NO ₃
M:	131.1 g/mol
CAS:	51-35-4
EINECS:	200-091-9
CS:	29339980

AppliChem GmbH

Ottoweg 4 • D-64291 Darmstadt • Phone +49 6151 9357 0 • Fax +49 6151 9357 11 • info.de@itwreagents.com • www.itwreagents.com
 CEO Joan Roget • Commerzbank Darmstadt • Bank 508 800 50 • Account 0186989900 IBAN DE24 5088 0050 0186 9899 00 • Swiftcode
 DRESDEFF508 • Finanzamt Darmstadt 007 225 00192 • Register court Darmstadt HRB Nr. 7340