

Specification

D-Cycloserine *BioChemica*

A1943

Physical Description:	Solid
Product Code:	A1943
Product Name:	D-Cycloserine <i>BioChemica</i>
Specifications:	Assay: min. 900 µg/mg $\alpha_{20^{\circ}\text{C/D}}; 5\%, 2\text{ N NaOH}$: +108° - +114° Solubility (10 %; H ₂ O): clear, colorless pH (10 %; H ₂ O; 20°C): 5.5 - 6.5 Residue on ignition: max. 0.5 % Loss on drying: max. 1 %
WGK:	1
Storage:	-20°C
Shipment:	RT
P Phrases:	P261 P281
Molecular Formula:	C ₃ H ₆ N ₂ O ₂
M:	102.10 g/mol
CAS:	68-41-7
EINECS:	200-688-4
CS:	29419000

AppliChem GmbH

Ottoweg 4 • D-64291 Darmstadt • Phone +49 6151 9357 0 • Fax +49 6151 9357 11 • info.de@itwreagents.com • www.itwreagents.com
 CEO Joan Roget • Commerzbank Darmstadt • Bank 508 800 50 • Account 0186989900 IBAN DE24 5088 0050 0186 9899 00 • Swiftcode DRESDEFF508 • Finanzamt Darmstadt 007 225 00192 • Register court Darmstadt HRB Nr. 7340