

Specification

Albumin (BSA) Fraction V (pH 5.2)

A2244

Physical Description:	Solid
Product Code:	A2244
Product Name:	Albumin (BSA) Fraction V (pH 5.2)
Specifications:	Assay (Protein): min. 96 % pH (10 %; H₂O; 20°C): 4.8 - 5.6 Protein (Albumin): min. 98 % NaCl: max. 5 % Water: max. 5 %
WGK:	1
Storage:	2 - 8°C
Origin:	from bovine serum
M:	approx. 68000 g/mol
CAS:	9048-46-8
EINECS:	232-936-2
CS:	35029070

AppliChem GmbH

Ottoweg 4 • D-64291 Darmstadt • Phone +49 6151 9357 0 • Fax +49 6151 9357 11 • info.de@itwreagents.com • www.itwreagents.com
 CEO Joan Roget • Commerzbank Darmstadt • Bank 508 800 50 • Account 0186989900 IBAN DE24 5088 0050 0186 9899 00 • Swiftcode DRESDEFF508 • Finanzamt Darmstadt 007 225 00192 • Register court Darmstadt HRB Nr. 7340