

Specification

Polyvinylpyrrolidone (K15) *BioChemica*

A2258

Physical Description:	Solid
Product Code:	A2258
Product Name:	Polyvinylpyrrolidone (K15) <i>BioChemica</i>
Specifications:	Assay (from N): 11.5 - 12.8 % Solubility (10 %; H₂O): clear, colorless to yellowish pH (5 %; H₂O): 3.0 - 5.0 Water (K.F.): max. 5 %
WGK:	1
Storage:	RT
M:	approx. 10000 g/mol
CAS:	9003-39-8
CS:	39059990

AppliChem GmbH

Ottoweg 4 • D-64291 Darmstadt • Phone +49 6151 9357 0 • Fax +49 6151 9357 11 • info.de@itwreagents.com • www.itwreagents.com
 CEO Joan Roget • Commerzbank Darmstadt • Bank 508 800 50 • Account 0186989900 IBAN DE24 5088 0050 0186 9899 00 • Swiftcode
 DRESDEFF508 • Finanzamt Darmstadt 007 225 00192 • Register court Darmstadt HRB Nr. 7340