

Specification

Tris for molecular biology

A2264

Solubility:	800 g/L (H ₂ O)
Physical Description:	Solid
Product Code:	A2264
Product Name:	Tris for molecular biology
Specifications:	DNases/RNases/Proteases: not detectable Assay (titr.): min. 99.9 % pH (1 M; H ₂ O; 20°C): 10.5 - 11.5 Water: max. 0.2 % Cu: max. 0.0001 % Fe: max. 0.0001 % Pb: max. 0.0001 % A (1 cm/1 M in H ₂ O) 430 nm: max. 0.004 260 nm: max. 0.03 280 nm: max. 0.02
WGK:	1
Storage:	RT
Molecular Formula:	C ₄ H ₁₁ NO ₃
M:	121.14 g/mol
CAS:	77-86-1
EINECS:	201-064-4
CS:	29221900

AppliChem GmbH

Ottoweg 4 • D-64291 Darmstadt • Phone +49 6151 9357 0 • Fax +49 6151 9357 11 • info.de@itwreagents.com • www.itwreagents.com
 CEO Joan Roget • Commerzbank Darmstadt • Bank 508 800 50 • Account 0186989900 IBAN DE24 5088 0050 0186 9899 00 • Swiftcode
 DRESDEFF508 • Finanzamt Darmstadt 007 225 00192 • Register court Darmstadt HRB Nr. 7340