

Specification

Heparin Sodium Salt *BioChemica*

A3004

Physical Description:	Solid
Product Code:	A3004
Product Name:	Heparin Sodium Salt <i>BioChemica</i>
Specifications:	Assay: min. 150 I.U. / mg Solubility (1 %; H₂O): clear, colorless pH (1 %; H₂O; 20°C): 5.5 - 8.0 Heavy metals (as Pb): max. 0.003 % Loss on drying: max. 8 % Na: 9.5 - 12.5 %
WGK:	1
Storage:	RT
Origin:	from porcine mucosa
M:	8000 - 25000 g/mol
CAS:	9041-08-1
CS:	30019091

AppliChem GmbH

Ottoweg 4 • D-64291 Darmstadt • Phone +49 6151 9357 0 • Fax +49 6151 9357 11 • info.de@itwreagents.com • www.itwreagents.com
 CEO Joan Roget • Commerzbank Darmstadt • Bank 508 800 50 • Account 0186989900 IBAN DE24 5088 0050 0186 9899 00 • Swiftcode
 DRESDEFF508 • Finanzamt Darmstadt 007 225 00192 • Register court Darmstadt HRB Nr. 7340