

Specification

L-Histidine base *BioChemica*

A3738

Solubility:	42 g/L (H ₂ O)
Physical Description:	Solid
Product Code:	A3738
Product Name:	L-Histidine base <i>BioChemica</i>
Headline Comment:	• All amino acids from AppliChem are of non-animal origin!
Specifications:	Assay (titr.): min. 99 % α20°C/D; 11 %, HCl 120 g/l: +11.8° - +12.8° Heavy metals (as Pb): max. 0.001 % Insoluble matter: passes test Other amino acids: max. 0.5 % Ammonium: max. 0.1 % Chloride: max. 0.02 % Sulfate: max. 0.02 %
WGK:	1
Storage:	RT
Molecular Formula:	C ₆ H ₉ N ₃ O ₂
M:	155.16 g/mol
CAS:	71-00-1
EINECS:	200-745-3
CS:	29332990

AppliChem GmbH

Ottoweg 4 • D-64291 Darmstadt • Phone +49 6151 9357 0 • Fax +49 6151 9357 11 • info.de@itwreagents.com • www.itwreagents.com
 CEO Joan Roget • Commerzbank Darmstadt • Bank 508 800 50 • Account 0186989900 IBAN DE24 5088 0050 0186 9899 00 • Swiftcode
 DRESDEFF508 • Finanzamt Darmstadt 007 225 00192 • Register court Darmstadt HRB Nr. 7340