

Specification

Phenol red sodium salt

A4233

Physical Description:	Solid
Product Code:	A4233
Product Name:	Phenol red sodium salt
Specifications:	Assay (photometr.): min. 95 % Loss on drying: max. 8 % Insoluble matter: max. 0.5 %
WGK:	1
Storage:	RT
Molecular Formula:	$C_{19}H_{13}NaO_5S$
M:	376.36 g/mol
CAS:	34487-61-1
EINECS:	252-057-8
CS:	29349990

AppliChem GmbH

Ottoweg 4 • D-64291 Darmstadt • Phone +49 6151 9357 0 • Fax +49 6151 9357 11 • info.de@itwreagents.com • www.itwreagents.com
 CEO Joan Roget • Commerzbank Darmstadt • Bank 508 800 50 • Account 0186989900 IBAN DE24 5088 0050 0186 9899 00 • Swiftcode
 DRESDEFF508 • Finanzamt Darmstadt 007 225 00192 • Register court Darmstadt HRB Nr. 7340