

Specification

Tris buffer pH 7.5 (1 M) for molecular biology

A4263

Physical Description:	Liquid
Product Code:	A4263
Product Name:	Tris buffer pH 7.5 (1 M) for molecular biology
Short Description:	additional product description: aqueous Solution
Specifications:	DNases/RNases/Proteases: not detectable pH (20°C; adjusted with HCl): 7.5 ± 0.2 Composition: Tris: 121.14 g/L (1 M)
WGK:	1
Storage:	RT
CS:	38221900

AppliChem GmbH

Ottoweg 4 • D-64291 Darmstadt • Phone +49 6151 9357 0 • Fax +49 6151 9357 11 • info.de@itwreagents.com • www.itwreagents.com
 CEO Joan Roget • Commerzbank Darmstadt • Bank 508 800 50 • Account 0186989900 IBAN DE24 5088 0050 0186 9899 00 • Swiftcode
 DRESDEFF508 • Finanzamt Darmstadt 007 225 00192 • Register court Darmstadt HRB Nr. 7340