


Specification

Pepsin

A4289


Physical Description:	Solid
Product Code:	A4289
Product Name:	Pepsin
Specifications:	Activity (Hemoglobin; pH 1.6; 25°C): min. 0.7 FIP-U/mg pH (2 %; H ₂ O): 3.5 - 5.0 Loss on drying (60°C; 4 h; 670 hPa; P ₂ O ₅): max. 5 %
Hazard pictograms	
WGK:	1
Storage:	2 - 8°C
Signal Word:	Danger
GHS Symbols:	GHS07 GHS08
H Phrases:	H315 H319 H334 H335
P Phrases:	P284 P305+P351+P338 P342+P311 P403+P233 P405 P501

AppliChem GmbH

Specification

Pepsin

A4289

Origin:	from porcine stomach
M:	~36000 g/mol
CAS:	9001-75-6
EINECS:	232-629-3
CS:	35079090
Index Nr.:	647-008-00-6
	

AppliChem GmbH