

Specification

Sodium Acetate anhydrous for molecular biology

A4555

Physical Description:	Solid
Product Code:	A4555
Product Name:	Sodium Acetate anhydrous for molecular biology
Specifications:	DNases/RNases/Proteases: not detectable Assay (titr.): min. 99 % pH (5 %; H ₂ O; 20°C): 7.5 - 9.5 Phosphate: max. 0.001 % Sulfate: max. 0.005 % Al: max. 0.001 % Ca: max. 0.001 % Cu: max. 0.0003 % Fe: max. 0.0005 % Mg: max. 0.0005 % Pb: max. 0.0005 % K: max. 0.02 % Zn: max. 0.0005 %
WGK:	1
Storage:	RT
Molecular Formula:	CH ₃ COONa
M:	82.03 g/mol
CAS:	127-09-3
EINECS:	204-823-8
CS:	29152900

AppliChem GmbH

Ottoweg 4 • D-64291 Darmstadt • Phone +49 6151 9357 0 • Fax +49 6151 9357 11 • info.de@itwreagents.com • www.itwreagents.com
 CEO Joan Roget • Commerzbank Darmstadt • Bank 508 800 50 • Account 0186989900 IBAN DE24 5088 0050 0186 9899 00 • Swiftcode
 DRESDEFF508 • Finanzamt Darmstadt 007 225 00192 • Register court Darmstadt HRB Nr. 7340