


Specification

Taq DNA Polymerase DNA-free

A5434

Product Code:	A5434
Product Name:	Taq DNA Polymerase DNA-free
Specifications:	<ul style="list-style-type: none"> : • Unit definition: see comment. : • supplied with reaction buffer (10X) Concentration: 5000 Units/ml Bacterial DNA: not detectable (min. 40 PCR cycle)
Storage:	-20°C
Shipment:	wet ice
Origin:	from <i>Thermus aquaticus</i> , recombinant (<i>E. coli</i>)
CS:	35079090
	

AppliChem GmbH

Ottoweg 4 • D-64291 Darmstadt • Phone +49 6151 9357 0 • Fax +49 6151 9357 11 • info.de@itwreagents.com • www.itwreagents.com
 CEO Joan Roget • Commerzbank Darmstadt • Bank 508 800 50 • Account 0186989900 IBAN DE24 5088 0050 0186 9899 00 • Swiftcode
 DRESDEFF508 • Finanzamt Darmstadt 007 225 00192 • Register court Darmstadt HRB Nr. 7340