


Specification

**Pepton Water, Buffered**

**A5841**

<b>Physical Description:</b>	Solid
<b>Product Code:</b>	A5841
<b>Product Name:</b>	Pepton Water, Buffered
<b>Specifications:</b>	<p>pH before autoclaving: approx. 7.2 (20°C)</p> <p><b>Composition:</b></p> <p>Sodium Chloride: 5.0 g/L</p> <p>Peptones: 10.0 g/L</p> <p>Buffers: 5.0 g/L</p>
<b>WGK:</b>	1
<b>Storage:</b>	RT
<b>CS:</b>	38221900
	

**AppliChem GmbH**

Ottoweg 4 • D-64291 Darmstadt • Phone +49 6151 9357 0 • Fax +49 6151 9357 11 • [info.de@itwreagents.com](mailto:info.de@itwreagents.com) • [www.itwreagents.com](http://www.itwreagents.com)  
 CEO Joan Roget • Commerzbank Darmstadt • Bank 508 800 50 • Account 0186989900 IBAN DE24 5088 0050 0186 9899 00 • Swiftcode  
 DRESDEFF508 • Finanzamt Darmstadt 007 225 00192 • Register court Darmstadt HRB Nr. 7340