

Specification

AppliCoat Plate Stabilizer

A7708

Product Code:	A7708
Product Name:	AppliCoat Plate Stabilizer
Specifications:	<ul style="list-style-type: none"> : • Solution for stabilizing coated surfaces in immunoassays : • Preservation of coated microtiter plates, polystyrene beads or glass slides : • preserves coated plates for 3 months up to 3 years depending on application : • free of proteins <p>pH: approx. 6.5</p>
WGK:	nwg
Storage:	2 - 8°C
H Phrases:	EUH208 EUH210
CS:	38221900

AppliChem GmbH

Ottoweg 4 • D-64291 Darmstadt • Phone +49 6151 9357 0 • Fax +49 6151 9357 11 • info.de@itwreagents.com • www.itwreagents.com
 CEO Joan Roget • Commerzbank Darmstadt • Bank 508 800 50 • Account 0186989900 IBAN DE24 5088 0050 0186 9899 00 • Swiftcode DRESDEFF508 • Finanzamt Darmstadt 007 225 00192 • Register court Darmstadt HRB Nr. 7340