


Specification

Protein Marker VI (10 - 245) prestained

A8889

Product Code:	A8889
Product Name:	Protein Marker VI (10 - 245) prestained
Specifications:	<ul style="list-style-type: none"> : • Three-color protein standard with 12 pre-stained. Proteins are covalently coupled with a blue chromophore except for two reference bands (one green and one red band at 25 kDa and 75 kDa respectively). : • Fragment size range in kDa: 10 - 245 : • Number of bands: 12
WGK:	2
Storage:	-20°C
Shipment:	wet ice
CS:	38221900
	

AppliChem GmbH

Ottoweg 4 • D-64291 Darmstadt • Phone +49 6151 9357 0 • Fax +49 6151 9357 11 • info.de@itwreagents.com • www.itwreagents.com
 CEO Joan Roget • Commerzbank Darmstadt • Bank 508 800 50 • Account 0186989900 IBAN DE24 5088 0050 0186 9899 00 • Swiftcode
 DRESDEFF508 • Finanzamt Darmstadt 007 225 00192 • Register court Darmstadt HRB Nr. 7340