

Specification

Penicillin - Streptomycin (100X) for cell culture

A8943

Product Code:	A8943
Product Name:	Penicillin - Streptomycin (100X) for cell culture
Short Description:	delivery form: stabilized, sterile filtrated, solution in 0.9 % NaCl
Specifications:	Endotoxins: passes test Composition: Penicillin: 10000 U/mL Streptomycin: 10 mg/mL
WGK:	1
Storage:	-20°C
Shipment:	wet ice in Germany, dry ice to abroad
H Phrases:	EUH208 EUH210
CS:	38221900

AppliChem GmbH

Ottoweg 4 • D-64291 Darmstadt • Phone +49 6151 9357 0 • Fax +49 6151 9357 11 • info.de@itwreagents.com • www.itwreagents.com
 CEO Joan Roget • Commerzbank Darmstadt • Bank 508 800 50 • Account 0186989900 IBAN DE24 5088 0050 0186 9899 00 • Swiftcode DRESDEFF508 • Finanzamt Darmstadt 007 225 00192 • Register court Darmstadt HRB Nr. 7340