

Specification

qPCR Mycoplasma Test Kit

A9019

Product Code:	A9019
Product Name:	qPCR Mycoplasma Test Kit
Specifications:	<ul style="list-style-type: none"> : ♦ lyophilized PCR Mix for the detection of mycoplasma in cell culture by real-time/qPCR : ♦ detects all of myoplasma species found in cell cultures : ♦ for 25 or 100 tests, respectively : ♦ Kit components <ul style="list-style-type: none"> : • PCR Primer Mix)incl. primers, probes, nucleotides) : • Positive template control : • Reaction Buffer Solution : • Water PCR grade : • internal control DNA
WGK:	nwg
Storage:	2 - 8°C
CS:	38221900

AppliChem GmbH

Ottoweg 4 • D-64291 Darmstadt • Phone +49 6151 9357 0 • Fax +49 6151 9357 11 • info.de@itwreagents.com • www.itwreagents.com
 CEO Joan Roget • Commerzbank Darmstadt • Bank 508 800 50 • Account 0186989900 IBAN DE24 5088 0050 0186 9899 00 • Swiftcode DRESDEFF508 • Finanzamt Darmstadt 007 225 00192 • Register court Darmstadt HRB Nr. 7340