

Specification

PBS tablets pH 7.4 (for 100 ml)

A9162

Product Code:	A9162
Product Name:	PBS tablets pH 7.4 (for 100 ml)
Specifications:	<ul style="list-style-type: none"> : • Ready to use in minutes : • Autoclaveable pH (25°C): 7.4 ± 0.05 <p>Composition:</p> <ul style="list-style-type: none"> KCL: 2.7 mM NaCl: 140 mM Phosphate: 10 mM
WGK:	1
Storage:	RT
Origin:	Tablets for preparation of 100 ml phosphate buffered saline (PBS)
CS:	38221900

AppliChem GmbH