

Specification

PBS tablets pH 7.4 (for 200 ml)

A9177

Product Code:	A9177
Product Name:	PBS tablets pH 7.4 (for 200 ml)
Specifications:	<ul style="list-style-type: none"> : • Ready to use in minutes : • Autoclaveable : • pH (25°C): 7.4 ± 0.05 <p>Composition:</p> <p>KCL: 2.7 mM</p> <p>NaCl: 140 mM</p> <p>Phosphate: 10 mM</p>
WGK:	1
Storage:	RT
Origin:	Tablets for preparation of 200 ml phosphate buffered saline (PBS)
CS:	38221900

AppliChem GmbH