


Specification

PCR Cyclor Validation Kit

A9742

Product Code:	A9742
Product Name:	PCR Cyclor Validation Kit
Short Description:	Usage: for conventional PCR block cyclors
Specifications:	two amplicons: Annealing temp. approx. 63°C one amplicon: Annealing temp. approx. 65°C no amplicon: Annealing temp. approx. 68°C
WGK:	nwg
Storage:	2 - 8°C
Shipment:	wet ice
CS:	38221900
	

AppliChem GmbH

Ottoweg 4 • D-64291 Darmstadt • Phone +49 6151 9357 0 • Fax +49 6151 9357 11 • info.de@itwreagents.com • www.itwreagents.com
CEO Joan Roget • Commerzbank Darmstadt • Bank 508 800 50 • Account 0186989900 IBAN DE24 5088 0050 0186 9899 00 • Swiftcode DRESDEFF508 • Finanzamt Darmstadt 007 225 00192 • Register court Darmstadt HRB Nr. 7340