

Specification

Tap (50 mm, PE) for 5 L PE canisters for sulfuric, hydrochloric and nitric concentrated acids. Also suitable for 10 L of Ethanol and Formaldehyde for clinical diagnosis **TP-05GRF**

Product Code:	TP-05GRF
Product Name:	Tap (50 mm, PE) for 5 L PE canisters for sulfuric, hydrochloric and nitric concentrated acids. Also suitable for 10 L of Ethanol and Formaldehyde for clinical diagnosis



AppliChem GmbH