



Biological Stains Special Promotion

In this promotion we present a selection of kits and dyes for different biological stains.

CE IVD Meet the requirements of the CE marking for in vitro diagnostic medical devices.

Kits ready to use

- All reagents are ready for use
- Supplied in easy-to-use dropper bottles
- Optimal staining of the sample
- Sufficient to perform up to 100 tests
- No additional equipment required



Product name	Application	Code	Package	Recommended list price	Promo recommended price
Gram-Hucker's Staining Kit for clinical diagnosis	The most widely used stain in microbiology to differentiate between Gram-positive and Gram-negative bacteria on the basis of their color retention	256649.0922	Kit	94.70 €	63.00 €
Kit for Fast Staining in Haematology (Fast Panoptic) for clinical diagnosis	Characterization of leukocytes. Results of a quality equal to classic staining methods (May Grünwald-Giemsa or Pappenheim) are obtained in only a few seconds	254807.0922	Kit	94.50 €	79.00 €
Masson's Trichrome Kit for clinical diagnosis	Indicated for connective tissue staining. It stains gametes, nuclei, nerve fibres, neuroglia, collagen, keratin and intracellular fibres. It can also be used to obtain a negative image of the Golgi apparatus	256692.0922	Kit Up to 100 tests	259.00 €	199.90 €
PAS Kit for clinical diagnosis	Periodic Acid-Schiff (PAS) is a staining method used to detect polysaccharides in formalin-fixed and paraffin embedded tissue sections	256676.0922	Kit Up to 100 tests	260.70 €	199.90 €
Reticulin Kit for clinical diagnosis	Reticulin is a mesh of fine fibers which provide support to the tissues. The Reticulin Kit is used for visualizing the presence of reticulin by impregnation with a silver salt	255115.0922	Kit Up to 50 tests	217.60 €	190.00 €

Powdered stains

For the preparation of solutions

Product name	Application	Code	Package	Recommended list price	Promo recommended price
Alcian Blue 8 GX (C.I. 74240) for clinical diagnosis	For histology, PAS-Alcian Blue staining, certified by the Biological Stain Commission	254584.1604	5 g	97.60 €	72.00 €
Eosin Yellowish (C.I. 45380) for clinical diagnosis	Vital staining and plasma staining	251299.1608	100 g	116.30 €	55.00 €
Hematoxylin 1-hydrate (C.I. 75290) for clinical diagnosis	Vaginal smear staining	251344.1604	5 g	82.50 €	65.00 €
		251344.1606	25 g	308.80 €	250.00 €



BO-027EN;201807

Stains in solution

Ready-to-use solutions for quick staining



Product name	Application	Code	Package	Recommended list price	Promo recommended price
Eosin Yellowish alcoholic solution 1% for clinical diagnosis	Solution for Hematoxylin-Eosin staining, in human and gynecological samples	256879.1210	500 ml	70.30 €	55.90 €
Giemsa's Azur-Eosin-Methylene Blue solution (slow) for clinical diagnosis	Blood smears and protozoa staining	251338.1610	500 ml	61.70 €	45.90 €
Harris Hematoxylin solution for clinical diagnosis	Solution for Hematoxylin-Eosin staining, in human and gynaecological samples	253949.1611	1000 ml	119.20 €	79.90 €
Harris Hematoxylin modified solution for clinical diagnosis	Solution for Hematoxylin-Eosin staining, in human and gynaecological samples, mercury free	256991.1610	500 ml	76.60 €	46.90 €
Kovacs' Reagent for clinical diagnosis	For detection of indole	252908.1608	100 ml	28.10 €	26.00 €
Lugol's Liquor with 0.33% of Iodine (diluted) for clinical diagnosis	For microscopy, bacteria staining according to Gram	256977.1609	250 ml	24.30 €	22.90 €
Lugol's Liquor with 0.4% of Iodine (diluted) for clinical diagnosis	For microscopy, bacteria staining according to Gram	251774.1608	100 ml	18.40 €	12.50 €
Lugol's Liquor with 5% of Iodine (concentrated) for clinical diagnosis	For microscopy, bacteria staining according to Gram	257041.1608	100 ml	24.30 €	12.75 €
May Grünwald's Eosin-Methylene Blue solution for clinical diagnosis	Blood smear staining	251416.1611	1000 ml	73.10 €	49.90 €
Papanicolaou's Solution EA 50 for clinical diagnosis	For cytology, cytoplasm staining	253594.1611	1000 ml	65.30 €	35.00 €
Papanicolaou's Solution OG 6 for clinical diagnosis	For cytology, cytoplasm staining of mature and keratinized cells	253892.1611	1000 ml	78.30 €	39.00 €



Products for microscopy

Products for different applications in microscopy such as mounting media, immersion oil and a product for paraffin cleaning.

Product name	Application	Code	Package	Recommended list price	Promo recommended price
DPX, mounting medium fast (toluene base) for clinical diagnosis	Non-aqueous mounting medium	255254.1610	500 ml	81.20 €	55.00 €
Eukitt®, mounting medium for clinical diagnosis	Adhesive and specimen preservative that can be used manually and in automated cover slipping equipment, fast drying	253681.0009	250 ml	82.50 €	55.00 €
Histofix® Spray fixative for clinical diagnosis	Fixative for Papanicolaou smears	256700.3408	6x100 ml	96.60 €	53.90 €
Immersion Oil for clinical diagnosis	Immersion oil for microscopy	251002.1208	100 ml	26.00 €	19.90 €
Paraffin Cleaner for clinical diagnosis	Microtomes cleaner used in the processing of human tissue	256876.3408	6x100 ml	171.90 €	118.00 €

Promotion valid from September 1 to December 31, 2018.

BO-027EN

AppliChem GmbH
 Ottoweg 4
 DE-64291 Darmstadt
 Germany
 Phone +49 6151 9357-0
 Fax +49 6151 93 57-11
info.de@itwreagents.com

Nova Chimica Srl
 Via G. Galilei, 47
 I-20092 Cinisello Balsamo
 (Milano) Italy
 Phone +39 02 66045392
 Fax +39 02 66045394
info.it@itwreagents.com

PanReac Química SLU
 Garraf, 2
 Polígono Pla de la Bruguera
 E-08211 Castellar del Vallès
 Barcelona (Spain)
 Phone +34 937 489 400
 Fax +34 937 489 401
info.es@itwreagents.com



www.itwreagents.com